

YOUNG[®]
SURGICAL

PuraExtend[®]
Resorbable
Collagen Membrane



Long-Term Barrier,
Lasting Support



Longer Resorption Time with Adequate Barrier Function

PuraExtend allows natural tissue coverage while acting as a barrier against epithelial growth. As healing progresses, soft tissue naturally covers the area within weeks. Some clinical experience has shown that PuraExtend provides an adequate barrier function when the membrane has been left exposed in the oral cavity.* Some individual results may vary.

With its **high mechanical strength and moderate confirmability**, the PuraExtend membrane can effectively address a host of clinical indications and surgical procedures[†], such as:

- Augmentation about implants placed in delayed and immediate extraction sockets
- Localized ridge augmentation for later implantation
- Alveolar ridge reconstruction for prosthetic treatment
- Filling of bone defects
- Guided bone regeneration in dehiscence defects
- Guided tissue regeneration procedures in periodontal defects



6–9 Month Resorption Time*

PuraExtend is a resorbable collagen membrane made from purified porcine dermis tissue, offering a balanced combination of strength, handling, and resorption. Its highly purified intact collagen provides **long-lasting protection for 6-9 months while maintaining barrier function.**



15mm x 20mm




20mm x 30mm



30mm x 40mm

*products
shown
at actual
size*



Dry Membrane Cross Section
50X Magnification

- Purified intact tissue matrix with high mechanical strength
- <0.01 micron pore size allows nutrients and growth factors while assuring cell barrier function
- Chemically crosslinked for longer in-vivo stability and minimizing potential immune response
- Significant tensile strength for ease of tacking and suturing

* Data on file at Regenity: CER.029-1-2-3

PuraExtend is Highly Conformable and Repositionable

Choose PuraExtend for long-lasting protection.

Item Number	Size
YSSPE15	15x20mm
YSSPE20	20x30mm
YSSPE30	30x40mm

To place your order, contact us at
youngsurgicalsolutions.com or call 1-877-540-3258



**ORDER
TODAY**





Supplies You Trust. Care That Counts.

© Young Surgical, 2025. All rights reserved.

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these names, logos, and brands does not imply endorsement.