

YOUNG[®]
SURGICAL

PuraFill[®]
Bone Graft
Composite Plug



Socket
Preservation,
Simplified



Similar to Natural Bone In Resorption and Remodeling Profile*


The PuraFill® Bone Plug is a simple, effective solution for filling extraction sockets and periodontal defects. Made from a synthetic calcium phosphate-based mineral with a carbonate apatite structure designed to be similar to natural bone. **Each plug consists of approximately 80% graft particulate combined with 20% Type I bovine collagen for enhanced handling and performance.**

Why Choose PuraFill:

- **Designed for Reliable Socket Preservation:**
Supports reliable socket preservation and consistent outcomes
- **Cost-Effective Solution:**
High-quality grafting at an affordable price.
- **Quick & Easy Application:**
Faster and simpler to use than bone particulate.
- **Hassle-Free Compliance:**
No donor tracking required.

Convenient Sizing and Packaging

PuraFill plugs are available in two standard sizes—**10 x 20mm** and **6 x 25mm**—ideal for various clinical needs. Each pack contains **10 plugs**, ensuring efficiency and convenience.




* Data on file at Regenity- In-vitro analysis

Expands to Fill Socket Site When Hydrated

PuraFill® can be placed pre-hydrated or dry, depending on clinician preference. It stabilizes within the grafting site, helping to contain the graft material and support clot retention for healing.

Once hydrated, PuraFill adapts to the shape of the defect, filling the socket for precise placement.



Scanning Electron
Micrograph of the
**PuraFill Plug at
50x magnification**

The Easy Way to Fill and Preserve

PuraFill is indicated for socket preservation, ridge augmentations and sinus lifts. It does not require a barrier membrane and its collagen component helps bind bone graft particles, reducing the risk of particle loss into the oral cavity.

With its ease of use and versatility, PuraFill provides a straightforward solution for socket and ridge preservation.



PuraFill Delivers Reliable, Adaptable Grafting

Choose PuraFill for simplified, cost-effective socket preservation.

Item Number	Size	Count
PFPLG1020	10x20mm	10 per box
PFPLG625	06x25mm	10 per box

To place your order, contact us at
youngsurgicalsolutions.com or call 1-877-540-3258



**ORDER
TODAY**





Supplies You Trust. Care That Counts.

© Young Surgical, 2025. All rights reserved.

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these names, logos, and brands does not imply endorsement.