


YOUNG[®]
SURGICAL

PuraFlex[®]

**Resorbable
Collagen Membrane**



**Highly Conformable
Tucks With Ease**




Easily Drapes Over the Defect Site for Superior Handling

PuraFlex® is a resorbable collagen dental membrane designed for strength and flexibility. Made from purified porcine peritoneum tissue, it is highly conformable and easy to handle.

Product Features*:

- Highly drapeable
- Excellent tensile strength
- Purified intact porcine peritoneum tissue
- GTR & GBR barrier functions



Scanning Electron
Micrograph of the
PuraFill Plug at
50x magnification

A Membrane That Bends, Not Battles

With high drapeability, PuraFlex adapts to the surgical site, allowing the corners to tuck smoothly under the gums for precise placement. Its balance of drapability and durability supports a seamless grafting process.

Why You'll Love PuraFlex:

- **Resorbable Barrier:**
Provides approximately 12–16 weeks of barrier function for guided tissue and bone regeneration (GTR & GBR).
- **Flexible Yet Strong:**
Minimal memory allows easy handling, while sufficient tensile strength supports suturing or tacking.
- **Highly Conformable:**
Adapts smoothly to the surgical site for precise placement.
- **Instrument-Friendly:**
Designed to resist sticking while maintaining stability for smooth handling during procedures




*Data on file at Regenity


The Right PuraGraft® Membrane for the Right Application

PuraGraft® membranes offer a range of conformability and resorption profiles to suit various clinical needs. Whether you need a long-lasting barrier or a highly drapable solution, we provide the right option for predictable and effective tissue regeneration.

PRODUCT CONFORMABILITY



PRODUCT RESORPTION RATES



Item Number	Size
YSSPF15	15x20mm
YSSPF20	20x30mm
YSSPF30	30x40mm

To place your order, contact us at
youngsurgicalsolutions.com or call 1-877-540-3258

PuraGraft Membranes Support the Following Dental Procedures:

- In augmentation around implants placed in immediate and delayed extraction sockets
- Localized ridge augmentation for later implantation
- Alveolar ridge reconstruction for prosthetic treatment
- Filling of bone defects
- Guided bone regeneration in dehiscence defects and guided tissue regeneration procedures in periodontal defects

With tailored resorption rates and varying levels of conformability, PuraGraft membranes ensure you have the right solution for every procedure.



**ORDER
TODAY**





Supplies You Trust. Care That Counts.

© Young Surgical, 2025. All rights reserved.

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these names, logos, and brands does not imply endorsement.